



# ***CMDR-ST Series***

## **Quick Setup Guide**





## Hardware Set-up

1. Uncoil main power cord from metal fingers on boom arm of CMDR transport stand. Unwinding of the power cord is required to avoid stress on the cord connection to the power supply box.




2. Raise tube head to desired position by depressing lever on boom arm with attention paid to the tension on the panel tether.



*Note: Consult owner's manual to rotate tube head for control panel to face laptop if preferred.*



2. Secure laptop computer to computer tray using the double-sided Velcro provided.
3. Plug in the CAT 5 cable that is attached to the computer power cable into the back of the laptop.
4. Plug the power adaptor cord in to the laptop.
5. Plug the power cord into grounded power source.
6. Turn on the computer.
7. Turn on the x-ray unit power switch located on top of the tube head.





## Software Initiation

1. The acquisition software will likely load several moments after the computer is powered up.

*If the acquisition software fails to load, double click the MinXray DX-R icon on the desktop to activate software.*



2. Login with User Name and Password (these will be set during installation).



Notes:

---

---

---



3. Input patient demographics by clicking the *new* button in the middle of screen. Click *insert patient* after entering information. Minimum fields to be populated are seen in **red**.

The screenshot shows the MINXRAY software interface. On the left, the 'patient' section contains input fields for title, first name (highlighted in yellow), last name, patient ID, DOB, age, sex (M/F), and study description. Below this is a section for study date/requesting physician with a date/time field showing '06/21/2018 15:42'. On the right, the 'Worklist' section has tabs for 'new', 'RIS', and 'emergency'. The 'new' tab is active, and a red arrow points to it. Below the tabs is a table with columns for last name, first name, study description, and accession number. The first row shows a patient with the time '15:42'. Below the table are 'insert patient' and 'cancel input' buttons.


*Note: Emergency Tab can be used when no patient demographics are supplied. This information can be added when available.*

This screenshot is similar to the previous one, but the 'emergency' tab in the Worklist section is selected, as indicated by a red arrow. The 'new' tab is no longer active. The rest of the interface, including the patient input fields on the left and the 'insert patient'/'cancel input' buttons, remains the same.






4. Select the *X-Ray* tab in upper right of screen to begin x-ray plan.



5. Select anatomy to be imaged on adult or child model on left of screen by clicking on the desired region to be imaged. Selected exams will populate the column on the right side of the screen. Use *edit* tab to delete unwanted or duplicate exams. Multiple exams can be selected before taking exposures.

*Note: MinXray DXR Imaging software contains imaging tutorials including correct positioning, sample digital radiographic images and video descriptions. To view the tutorials, click on the position thumbnail of selected study.*






6. Once all exams are planned, select *exposure* tab.



7. Select the position to be imaged, it will be highlighted when selected.





## X-Ray Exposure

1. Remove imaging panel from storage compartment located on the frame of the CMDR Stand.



2. Uncoil panel tether cord from storage hooks and extend to desired length.  
Be aware of cable placement when positioning to avoid damage to panel or cable.
3. Position imaging panel with logo side up or facing patient.







4. Adjust focal spot of tube head 40 inches (100cm) from, and perpendicular to panel using tape measure located on the bottom of collimator.

*Note: For extra tube head height, see owner's manual for frame extension.*

5. Turn on collimator light using button on front of collimator and narrow the light field to desired imaging field. Collimator light will automatically turn off after 25 seconds.



6. Set appropriate exposure values according to body part and exposure chart located on the back of tube head using kV and mAs/seconds touchpad on front of x-ray tubehead.

***Note: The LED readout for mAs/seconds toggles between mAs and seconds as indicated by the indicator light under each value. The display can be changed between mAs and seconds by pressing the # button and the down arrow below the kV display.***



7. Using exposure switch, press and hold button until exposure light and audible indication is complete.



The exposure switch has a two-stage process. The first stage pre-heats the filament of the x-ray tube. An x-ray exposure is made when the second stage is fully depressed.

Press the exposure button half way down until the green ready light is illuminated on the front of the tube head, then press button all the way down until orange X-ray light illuminates for x-ray exposure. Hold the button down until the audible tone stops, or the orange light goes out, then release button.


**Premature release of the exposure button may result in an incomplete exposure.**

## Post-exposure Image Processing


1. Manipulate image for proper display and measurements using QA tools under image.




2. Selecting the *Eyeball* will bring up a new window with advanced tools and annotations which will be covered in training.




3. Delete image (select *trash can*) for improper position, artifacts or motion and retake.



4. Select *check mark* under image to accept image and move to next position.




5. Once all images are taken select *finish study* check mark.






## Viewing Complete Studies

1. To review studies, select the Lists tab on upper right screen and select study to be viewed.




2. Studies highlighted in green are complete, studies in red are incomplete or have not been QA-ed or finished.
3. Select *Create CD* on bottom right of screen to burn images to CD.






## Shut Down / Storage Procedures

1. Click on Operator key in upper left of screen.




2. Click on Exit tab in lower right to exit DX-R software.





3. Click on start and then select Turn Off Computer.
4. Turn off X-Ray unit using power switch located on top of x-ray tube.
5. Place DR panel in storage compartment and carefully wrap panel cord around storage hooks.
6. Close and latch panel storage compartment door. Be sure that the panel cable is correctly placed and is not pinched in storage compartment door.



7. Lower x-ray tube by depressing gas spring actuator on boom arm of CMDR stand.
8. Lower unit until tube head gently rests against rubber pad on frame.
9. Wrap power cord around storage hooks on boom arm of CMDR stand.
10. If laptop is to be removed, unplug power adaptor and panel connector before removing from holder tray.
11. Fold down holder tray and secure in down position with latch.

Technical Service is available by phone: 800-221-2245 or 847-564-0323  
E-mail: [drimaging@minxray.com](mailto:drimaging@minxray.com) during the hours of 8:00 a.m. to 4:30 p.m. Central Time,  
Monday through Friday, except holidays.  
Telephone and online support 24/7 support is available as a Service Contract option.